


Wondfo



CE CFDA


Optical Coagulation Analyzer

Optical Coagulation Analyzer OCG-102


The first POCT coagulation analyzer that can test
PT/INR, APTT, TT, FIB and ACT

Portable	Small size (weight only 850g); Use with rechargeable battery	Rapid	Get results within 3 mins ; Short turn around time
Convenient	One step operation, no accessory required.	Accurate	CV = 5%
Sample Type	Citrated venous whole blood, capillary whole blood	Sample Volume	20μL


SIMPLE TEST PROCEDURE



Inserting



Dropping



Reading

- Disposable Dried-reagent test strips
- Room temperature storage
- Long shelf-life: 18 months
- PT/INR can be performed with fresh fingertip capillary whole blood
- One drop of blood for monitoring oral anticoagulation therapy (warfarin)



- Screening test to assess the exogenous and endogenous coagulation system function
- Monitoring anticoagulant treatment
- Monitoring thrombolytic therapy
- Monitoring the patient with certain disease, such as hepatopathy, infection, burn injury

